

SECTION 1 — SAFETY AND INTENDED PURPOSE

The SimpStrap Pro is a Class I passive wound dressing intended for management of minor skin wounds in community settings. The device has been designed to ensure that it does not compromise the health and safety of patients, users, or third parties when used as intended. A risk management process has been conducted in accordance with ISO 14971:2019. All identified risks have been reduced to acceptable levels.

SECTION 2 — DESIGN AND CONSTRUCTION

The device consists of woven cotton gauze and medical grade adhesive. All materials have been selected for their established history of use in wound care. Biocompatibility has been assessed with reference to ISO 10993-1:2018. No materials of biological origin are used. The device is not flammable under normal conditions of use.

SECTION 3 — INFORMATION PROVIDED WITH THE DEVICE

The device is supplied with labelling compliant with the requirements of Schedule 1 of the Therapeutic Goods (Medical Devices) Regulations 2002. The label includes:

- Device name and intended purpose
- Manufacturer name and address
- Batch/lot number and expiry date
- Instructions for use
- Single use only indication
- Storage conditions

All information is provided in English and uses the standard ISO 15223-1 symbols.

SECTION 4 — MANUFACTURING AND QUALITY

Manufacturing is conducted under ISO 13485:2016 certified quality management system.

Shelf life has been validated to 36 months under specified storage conditions.

Transport and storage studies demonstrate the device is not adversely affected by transport or storage conditions within the validated range.

SECTION 5 — BENEFIT-RISK ANALYSIS

Clinical benefits of wound management with SimpStrap Pro outweigh any residual undesirable effects. No serious adverse events are anticipated for a Class I passive wound dressing used as intended.